



AimStrip® Hb

HEMOGLOBIN (Hb) TEST SYSTEM

**CLIA
WAIVED**



RAPID HEMOGLOBIN & HEMATOCRIT RESULTS!

For the quantitative determination of hemoglobin in non-anticoagulated capillary whole blood or anticoagulated venous whole blood. Estimation of hematocrit is provided for normal hemoglobin values.¹



"AimStrip® Hemoglobin provides lab quality accuracy for Hb and Hct testing!"



LOOK AT THESE GREAT FEATURES!

One Test. Two results!

1 HEMOGLOBIN
Measurement range
5.6 - 23.5g/dL

2 HEMATOCRIT
Estimated for normal
hemoglobin values.¹

Rapid

- ⇒ Simply insert test, apply specimen with capillary tube & read results <15 sec.
- ⇒ No specimen preparation required

POC Testing

- ⇒ Handheld meter for POC testing
- ⇒ Easy to read LCD results

Minimal Training Time

- ⇒ User friendly meter makes testing simple

Zero Maintenance

- ⇒ No servicing of meter is required
- ⇒ Meter has 3 year warranty

Convenient

- ⇒ No specimen preparation required
- ⇒ Uses batteries or optional AC adaptor

CAP / CLSI Procedures

- ⇒ Already written, available as MS Word & pdf.
- ⇒ Just add your information.

CPT CODE: 85018 QW.....Reimbursement \$3.31*



www.germainelabs.com



800.854.8446
210.692.4192
210.692.4198 Fax



**THE TESTS YOU WANT.
THE ACCURACY YOU NEED!**

DISTRIBUTED BY:

¹ See Package Insert for complete procedure.
*CPT Code is assigned by Medicare. Check with each state for exact reimbursement.

AimStrip® Hb

HEMOGLOBIN (Hb) TEST SYSTEM

**One Test.
Two results!**



PRODUCT NAME

AimStrip® Hb Starter Kit

REF#

78200

- ⇒ 1 Meter
- ⇒ 200 Test Cartridges
- ⇒ 240 Capillary Tubes
- ⇒ 2 Optical Verifiers
- ⇒ 4 Code Chips
- ⇒ 1 Carrying Case
- ⇒ 4 AAA Batteries
- ⇒ 1 User Manual
- ⇒ 4 Results Record Pads

AimStrip® Hb Meter

78101

- ⇒ 1 Meter
- ⇒ 2 Optical Verifiers
- ⇒ 1 Carrying Case
- ⇒ 4 AAA Batteries

AimStrip® Hb Test Combo

78250

- ⇒ 50 Test Cartridges
- ⇒ 60 Capillary Tubes

AimStrip® Hb Optical Verifiers

78102

- ⇒ 2 Optical Verifiers

AimStrip® Hb Control Set

78103

- ⇒ 1x2mL Control Solution 0
- ⇒ 1x2mL Control Solution 1
- ⇒ 1x2mL Control Solution 2

AimStrip® Hb AC Adaptor

78111

- ⇒ 1 AC Adaptor for Meter



CLIA
WAIVED



PRODUCT SPECIFICATIONS

Methodology	Reflectance Photometry
Calibration	Automatic & Factory Calibrated
Sample Type	Capillary and Venous Whole Blood
Sample Size	10µL
Hb Measurement Range	5.6 to 23.5 g/dL
Hematocrit	Estimated, see package insert
Test Time	<15 Seconds
Battery Life	2,700 Tests or 360 Hours
Power Source	4 AAA Batteries or AC Adaptor
Automatic Shutoff	8 minutes
Memory	1,000 Tests, Date/time & ID Number

Bottom Line
= Accurate, Rapid, Easy!

THE TESTS YOU WANT. THE ACCURACY YOU NEED!

GERMAINE
LABORATORIES, INC.

www.germainelabs.com

