



AIMSTRIP® URINE ANALYZER 2

**CLIA
WAIVED**



*AimStrip® Urine Analyzer 2
REF #51222*

AimStrip® Urine Analyzer 2 is for the semi-quantitative detection of the following analytes in urine: Glucose, Bilirubin, Ketone (Acetoacetic acid), Specific Gravity, Blood, pH, Protein, Urobilinogen, and Leukocytes as well as the qualitative detection of Nitrite.

**Urinalysis
Made Easy!**

LOOK AT THESE GREAT FEATURES!

- ⇒ Single Test Mode, reads 1 strip at a time¹
- ⇒ Continuous Test Mode, reads 120 strips
- ⇒ Automatic calibration for accurate results
- ⇒ Automatically saves the last 2000 results
- ⇒ Reads 2, 3, 4, 5, 6, 7, 8, 9 and 10 parameter strips
- ⇒ Includes RS232C and USB ports to export results

PROGRAM LANGUAGES: English, Español, Portuguese, Français, Italiano.



PAPER	SIZE	REF#
Printer Paper	4 rolls/pk	51224

- ⇒ Thermal paper to print 200 results per roll

Printer Peel & Stick Paper	4 rolls/pk	51234
---------------------------------------	------------	-------

- ⇒ Peel and stick printer paper to print 100 results per roll
- ⇒ Conveniently place results into patient records



**STORE AT
2° TO 8°C**

CONTROLS

UA Liquid Controls	SIZE	REF#
	2 x 5mL	51115

- ⇒ Negative control (1 x 5mL), Positive control (1 x 5mL)

UA Liquid Controls	6 x 10mL	51161
---------------------------	----------	-------

- ⇒ Negative control (3 x 10mL), Positive control (3 x 10mL)

UA Liquid Controls Plus	2 x 5mL	51116
--------------------------------	---------	-------


- ⇒ Negative control (1 x 5mL), Positive control (1 x 5mL)
- ⇒ For pregnancy tests, AimTab™ Bilirubin Tablets, AimTab™ Ketone Tablets, and AimTab™ Reducing Substances Tablets & UA strips

UA DipTube Controls	2 x 12mL	51112
----------------------------	----------	-------

- ⇒ Negative control (1 x 12mL), Positive control (1 x 12mL)

UA DipeTube Controls	4 x 12 mL	51142
-----------------------------	-----------	-------

- ⇒ Negative control (2 x 12mL), Positive control (2 x 12mL)



OPTIONAL PRODUCT

- Barcode Reader #51320**
- ⇒ Optional to record patient ID



www.germainelabs.com



800.854.8446
210.692.4192
210.692.4198 Fax

**THE TESTS YOU WANT.
THE ACCURACY YOU NEED!**

DISTRIBUTED BY:

