

AimTab™ Bluing Tablets

AimTab™ Bluing Tablets are for use in all drug testing collections to prevent adulteration and for leak detection by plumbers



Quickly Turns Water Blue!



For a Deep Vibrant Blue

- AimTab™ Bluing Tablets, 100/bottle, Ref #40021
- Place tablet into toilet bowl before collection and turns water vibrant blue in less than a minute
- Use one tablet per flush
- Meets Federal DOT Regulation 49 CFR 40.41
- Plumbers use for testing customers toilets for leaking flapper and ball valves by placing into toilet tank
- CAP / CLSI formatted procedures available via fax or email - just add your name!
- Store at 10° to 35°C (50° to 95°F)

Related Products:

- Aim™ Bluing Powder, 4 oz. Ref #40022
- Aim™ Bluing Liquid, 4 oz. Ref #40023

Tests You Want,
Accuracy You Need!



www.germainelabs.com

800.854.8446

210.692.4192

210.692.4198 Fax

*See Package Insert for complete procedure. Rev. 7-16
©2016 Germaine Laboratories, Inc., CSM 0800-59

Distributed by: