

AimTab™ BILIRUBIN TABLETS



CLIA
WAIVED

AimTab™ Bilirubin Tablets are for the confirmatory detection of bilirubin in urine. This test is for in-vitro diagnostic use.

**Bilirubin Confirmation
Made Easy!**

Same test procedure as Ictotest®, a registered trademark of Siemens Healthcare.

LOOK AT THESE GREAT FEATURES!

AimTab™ Bilirubin Tablets,
REF #13200


100 TABLETS
PER BOTTLE



READING RESULTS¹

**POSITIVE RESULTS =
BLUE OR PURPLE COLOR**


The test is **POSITIVE** IF A **BLUE OR PURPLE COLOR DEVELOPS** ON THE TEST MAT. THIS INDICATES THAT **BILIRUBIN IS PRESENT**. A **slight pink or red color** should be ignored.



**NEGATIVE RESULTS =
NO BLUE OR PURPLE COLOR**

The test is **NEGATIVE** IF **NO BLUE OR PURPLE COLOR DEVELOPS** ON THE TEST MAT.



 CAP / CLSI formatted procedures available via fax or email - just add your name!



www.germainelabs.com

GERMAINE[®]
LABORATORIES, INC.



800.854.8446
210.692.4192



210.692.4198 Fax



Provides quick & easy urine bilirubin confirmatory results in **60 seconds**.



Detects 0.05 to 1.0 mg/dL
(0.8 to 1.7 μ mol/L).



Confirms "urine strip" results.



Eliminates the need to perform blood test for bilirubin.



Store at room temperature.

**THE TESTS YOU WANT.
THE ACCURACY YOU NEED!**

DISTRIBUTED BY:

