



# AimTab™ KETONE TABLETS

**CLIA  
WAIVED  
FOR  
URINE**



AimTab™ Ketone Tablets, REF #13100

**100 TABLETS  
PER BOTTLE**



## RELATED PRODUCTS

### CONTROLS

**Ketone Serum Controls**      **SIZE** 2x2mL      **REF#** 13112

- ⇒ Negative control (1x2mL)
- ⇒ Positive control (1x2mL)
- ⇒ For use with Ketone tablet tests
- ⇒ Store at 2° to 8°C



**UA Liquid Controls**      **SIZE** 2x5mL      **REF#** 51116

- ⇒ Level 1 (1x5mL)
- ⇒ Level 2 (1x5mL)
- ⇒ Used for Ketone control
- ⇒ Store at 2° to 8°C



[www.germainelabs.com](http://www.germainelabs.com)

**GERMAINE**  
LABORATORIES, INC.



800.854.8446

210.692.4192





210.692.4198 Fax

AimTab™ Ketone Tablets are for the semi-quantitative determination of ketones (acetoacetic acid and acetone) in urine, serum, plasma and whole blood.

**Ketone Testing  
Made Easy!**

Same test procedure as Acetest®, a registered trademark used by Bayer HealthCare LLC.

## READING RESULTS¹

POSITIVE RESULTS = PURPLE COLOR		NEGATIVE RESULTS = NO PURPLE COLOR	
	The test is POSITIVE IF A <b>PURPLE COLOR DEVELOPS</b> ON THE TABLET	The test is NEGATIVE IF <b>NO PURPLE COLOR DEVELOPS</b> ON THE TABLET	
<b>READ</b>			
urine at 30 seconds		serum & plasma at 2 minutes	whole blood at 10 minutes



## Look at these great features:

- ⇒ Quick & easy ketone results for testing urine, serum, and whole blood.
- ⇒ Detects as little as 5 mg/dL of acetoacetic acid in urine.
- ⇒ Specific for acetoacetic acid and acetone and will not react with betahydroxybutyric acid.
- ⇒ CLIA moderately complex when testing serum, plasma or whole blood.
- ⇒ CAP / CLSI formatted procedures available via fax or email - just add your name!
- ⇒ Store at room temperature.

**THE TESTS YOU WANT.  
THE ACCURACY YOU NEED!**

**DISTRIBUTED BY:**



¹ See Package Insert for complete procedure.