



# Compliance<sup>®</sup>Gold<sup>®</sup> iFOB Test

A rapid chromatographic assay for the qualitative detection of fecal occult blood in human stool samples.

**Compliance  
Saves Lives!™**

**Compliance Gold® iFOB makes testing for occult blood more convenient for you and easier for your patients.**

**CLIA  
WAIVED**



## \*\* COMPLIANCE MONITORING SYSTEM



### STEP 1

Patient reads HIPAA Compliant Patient Agreement Form, signs form and enters their phone number.



### STEP 2

Doctor/physician faxes the Patient Agreement Form to the secure fax line at 888-797-4362.



### STEP 3

A series of reminder messages will remind patients to return completed tests to their doctor.



[www.germainelabs.com](http://www.germainelabs.com)

**GERMAINE<sup>®</sup>**  
LABORATORIES, INC.



800.854.8446

210.692.4192



210.692.4198 Fax

©2019 Germaine Laboratories, Inc., CSM 0800-02, Rev. 04-19



## COMPLIANCE MONITORING\*\* AT NO COST

Our team provides a HIPAA Compliant program that reminds patients to return Collection Kits.

### Look at these great features:

- ➡ **Postage Provided for Return Mailers**  
Patient Friendly! No more trips to the post office to buy postage.
- ➡ **Return Labels Included**  
Ensures tests are returned to the proper address.
- ➡ **Fits in Mail slot**  
Easier for patients to RETURN from the convenience of their own home!
- ➡ **Patient Instructions online 24/7**  
Always available for patients to reference for the collection and return process.
- ➡ **All Components Pre-Packaged in Mailer**  
Everything your Patients Need in one Return Mailer

CPT CODE: 82274QW - G0328QW.....Reimbursement \$22.22\*

**THE TESTS YOU WANT.  
THE ACCURACY YOU NEED!**

**DISTRIBUTED BY:**



\* See Package Insert for complete procedure.

\* CPT codes are assigned by Medicare for reimbursement.


Check online at <http://www.cms.hhs.gov> or ask us for your state's reimbursement.

# ComplianceGold® iFOB Test

IMMUNOLOGICAL FECAL OCCULT BLOOD TEST

**Compliance  
Saves Lives!™**

## SIMPLE PROCEDURE¹



CLIA  
WAIVED



30 TESTS PER KIT

**DETECTS:**  
Human hemoglobin In feces.

PRODUCT	REF#
Compliance Gold® iFOB Test	74034
Compliance Gold® iFOB Patient Collection Kits	74014
15 Complete Return Mailers	

READ RESULTS IN 5 MINUTES.

## READING RESULTS¹

VALID TEST RESULTS			INVALID TEST RESULTS	
NEGATIVE	POSITIVE	POSITIVE	INVALID	INVALID

FEATURE	Compliance Gold®	Guaiaac
Dietary Restrictions	No restrictions	Dietary Restrictions
Medication Restrictions	48 hrs some medications, none on Vitamin C	7 days some medications, 3 days Vitamin C over 250mg
Number of Samples	1 required	Up to 3 required
Reimbursement	\$22.22*	\$4.54*
Sensitivity	0.05 µg/ml of human hemoglobin in feces	90 µHb/ml of hemoglobin in feces
Specificity	Specific to human hemoglobin	Detects peroxidase activity and animal hemoglobin

**BOTTOM  
LINE .....**



IMPROVED  
COMPLIANCE!



Behind the curve  
of medical progress.

## ECONOMY VERSION

The following kits do not include Compliance Monitoring and Postage.

PRODUCT	REF#
Compliance Gold® iFOB Simple	74030
Compliance Gold® iFOB Simple Patient Collection Kits	74015
15 Complete Return Mailers	

## BASIC VERSION

The following version includes test device & collection bottle only.

PRODUCT	REF#
Compliance Gold® iFOB Easy	74033

## CONTROLS

Positive & Negative controls.

PRODUCT	SIZE	REF#
iFOB Controls	2 x 1mL	74211

¹ See Package Insert for complete procedure.

\* CPT codes are assigned by Medicare for reimbursement.

Check online at <http://www.cms.hhs.gov> or ask us for your state's reimbursement.

**Our fecal occult blood products** are manufactured to the highest quality standards utilizing current Good Manufacturing Practices (cGMP) as overseen by the FDA and the International Standardization Organization (ISO) 13485:2016. These certifications attest to the quality standards set and achieved for the manufacture of our fecal occult blood products.

Certified training is available for all registered visitors at: [www.compliancegold.org](http://www.compliancegold.org).

**THE TESTS YOU WANT. THE ACCURACY YOU NEED!**



[www.germainelabs.com](http://www.germainelabs.com)

