



FEBRILE PRODUCTS

Febrile Antigens are intended for the detection of bacterial agglutinins in serum associated with a bacterial infection.



**STORE AT
2° TO 8°C**

The slide test is recommended as a screening procedure to establish the presence or absence of homologous antibodies. The tube test is intended to establish antibody titer. The Febrile Controls are intended for use as positive & negative controls for these tests.

For testing patients serum using slide or tube test procedures



CONTROLS

	SIZE	REF#
Brucella Abortus Antiserum <i>Brucella control</i>	3mL	1100
Febrile Positive Control Positive control for <i>B. abortus</i> , <i>P. Ox19</i> , <i>Sal. H a</i> , <i>Sal. H b</i> , <i>Sal. H d</i> , <i>Sal. O Gr. D</i>	5mL	1600
Febrile Negative Control <i>For all febrile antigens</i>	5mL	1601
Francisella Tularensis Antiserum <i>F. tularensis control</i>	2mL	1702
Proteus Poly Antiserum Positive control for <i>Proteus OX2</i> , <i>OX19</i> , <i>OXK</i>	3mL	1303
Salmonella H Poly Positive Control Positive control for <i>Sal. H a,b,c,d</i>	3mL	1606
Salmonella O Poly Positive Control Positive control <i>Sal. O Gr. A, B, C, D</i>	3mL	1611

BACTERIAL ANTIGENS

	SIZE	REF#	
Brucella Abortus Antigen	5mL	1000	Slide & tube test
Brucella Melitensis Antigen	5mL	1002	Slide & tube test
Francisella Tularensis Antigen	3mL	1703	Slide & tube test
Proteus OX2 Antigen	5mL	1201	Slide & tube test
Proteus OX19 Antigen	5mL	1202	Slide & tube test
Proteus OXK Antigen	5mL	1200	Slide & tube test
Salmonella H a Antigen (paratyphoid A)	5mL	1400	Slide & tube test
Salmonella H b Antigen (paratyphoid B)	5mL	1401	Slide & tube test
Salmonella H c Antigen (paratyphoid C)	5mL	1402	Slide & tube test
Salmonella H d Antigen (tuphoid H)	5mL	1403	Slide & tube test
Salmonella O Gr. A Antigen (somatic ag.)	5mL	1500	Slide & tube test
Salmonella O Gr. B Antigen (somatic ag.)	5mL	1501	Slide & tube test
Salmonella O Gr.C Antigen (somatic ag.)	5mL	1502	Slide & tube test
Salmonella O Gr. D Antigen (typhoid O)	5mL	1503	Slide & tube test



⇒ CAP / CLSI formatted procedures available via fax or email - just add your name!



www.germainelabs.com



800.854.8446
210.692.4192
210.692.4198 Fax



**THE TESTS YOU WANT.
THE ACCURACY YOU NEED!**

DISTRIBUTED BY: