

PREGNANCY TESTS

THE FORMAT YOU WANT; THE SENSITIVITY YOU NEED!



**CLIA
WAIVED**
FOR
URINE

CAP / CLSI Procedures

- ⇒ Available via email or fax
- ⇒ Just add your name and approvals

CONTROLS

HCG Urine Controls **SIZE** 2x5mL **REF#** 71002

For use with urine pregnancy tests

- ⇒ Negative control (1 x 5mL)
- ⇒ Positive Control (1 x 5mL)
- ⇒ Store at 2° to 8°C



HCG Serum Controls **SIZE** 3x2.5mL **REF#** 71003

For use with serum pregnancy tests

- ⇒ Negative control (1 x 2.5mL)
- ⇒ Low Positive Control (1 x 2.5mL)
- ⇒ High Positive Control (1 x 2.5mL)
- ⇒ Store at 2° to 8°C



Pregnancy Testing Made Easy!

AimStep® PREGNANCY

- Urine cassette method
- Sensitive to 20 mIU/mL hCG
- Results in 3 minutes
- Store at room temperature
- CLIA Waived, also available in OTC

AimStep® COMBO PREGNANCY

- Serum or Urine cassette method
- Sensitive to 20 mIU/mL hCG
- Urine results in 3 minutes
Serum results in 5 minutes
- Store at room temperature
- CLIA Waived for urine testing

SIZE	30 TESTS	30 TESTS
REF#	97030	97730



AimStrip® PREGNANCY

- Urine dipstick method
- Sensitive to 20 mIU/mL hCG
- Results in 3 minutes
- Store at room temperature
- CLIA Waived, also available in OTC

SIZE	30 TESTS
REF#	88530

CPT CODE URINE: 81025 QW.....Reimbursement \$8.84*
CPT CODE SERUM: 84703Reimbursement \$10.49*



www.germainelabs.com



800.854.8446
210.692.4192
210.692.4198 Fax



THE TESTS YOU WANT.
THE ACCURACY YOU NEED!

DISTRIBUTED BY: