



SafeCup™ III

SafeCup™ III tests are rapid qualitative tests for the detection of drugs of abuse in urine.

SafeCup™ provides certified accurate results that are easy to read and require minimal training.

Quality Performance in your price range!

CLIA
WAIVED

25 CUPS
PER CASE



Store at 2° to 30°C

DRUG CUT-OFF LEVELS

Amphetamine	AMP	1000 ng/mL
Benzodiazepine	BZO	300 ng/mL
Barbiturates	BAR	300 ng/mL
Buprenorphine	BUP	10 ng/mL
Cocaine	COC	300 ng/mL
Ecstasy	MDMA	500 ng/mL
Marijuana	THC	50 ng/mL
Methadone	MTD	300 ng/mL
Methamphetamine	MET	1000 ng/mL
Morphine	MOP	300 ng/mL
Opiates	OPI	2000 ng/mL
Oxycodone	OXY	100 ng/mL
Phencyclidine	PCP	25 ng/mL
Propoxyphene	PPX	300 ng/mL
Tricyclic Antidepressants	TCA	1000 ng/mL

* OTHER CUT-OFF LEVELS ARE AVAILABLE
UPON REQUEST



www.germainelabs.com



800.854.8446
210.692.4192
210.692.4198 Fax

©2019 Germaine Laboratories, Inc., CSM-0800-04, Rev 03-19

LOOK AT THESE GREAT FEATURES!



The Total Collection, Testing,
and Transport System



Easier to Use, Easier to Read



Adheres to the SAMHSA
Cutoff Levels



Includes "Click seal" in cap
for improved leak resistance



Temperature Strip Included



Quicker Results - Read
Results in 5 minutes

**THE TESTS YOU WANT.
THE ACCURACY YOU NEED!**

DISTRIBUTED BY:



SafeCup™ III ALL-IN-ONE DRUG CUP

25 CUPS PER CASE




SIMPLE PROCEDURE




READING RESULTS

EXAMPLE OF VALID RESULTS
5 PANEL CUP



EXAMPLE OF VALID RESULTS
10 PANEL CUP



CLIA
WAIVED

STORE AT
2° TO 30°C

PRODUCT	REF#	DETECTS
MultiDrug 5	45475	COC, AMP, MET, THC, OPI
MultiDrug N5	45485	COC, THC, OPI, AMP, PCP
MultiDrug 10	40105	COC, THC, OPI, AMP, MET, PCP, BZO, BAR, BAR, MTD, MDMA
MultiDrug 12	41513	AMP, BAR, BZO, COC, THC, MDMA, MTD, MET, OPI, OXY, PCP, BUP
MultiDrug 12	41577	AMP, BAR, BUP, BZO, COC, MET, MDMA, MOP, MTD, OXY, PCP, THC
MultiDrug 12	41578	w/3AD (OX/SG/PH) COC, AMP, MET, THC, MOP, MTD, MDMA, OXY, PCP, BAR, BZO, BUP
MultiDrug 14	44418	AMP, BAR, BUP, BZO, COC, MET, MDMA, MOP, MTD, OXY, PCP, PPX, TCA, THC

OTHER CONFIGURATIONS ARE AVAILABLE UPON REQUEST

CONTROLS

DOA Positive Control

SIZE 1 x 5mL REF# 40002

Positive Control
(AMP, BAR, BUP, BZO, COC, MDMA, MET, MOP, MTD, OPI, OXY, PCP, PPX, TCA, THC)

⇒ Store at 2° to 8°C



DOA Negative Control

SIZE 1 x 5mL REF# 40003

Negative control
for all analytes

⇒ Store at 2° to 8°C




AimStrip® Adulteration 7

SIZE 25 tests REF# 40007

Urine dipsticks for
Creatinine, Nitrite, Glutaraldehyde, pH, Specific Gravity, Oxidant, PCC

⇒ Store at room temperature



Germaine Laboratories, Inc. has been at the leading edge of on-site instant testing technology since 1990. Many features from our other product developments have been incorporated into our new SafeCup™ III Drugs of Abuse testing products. Our Drugs of Abuse products are manufactured to the highest quality standards utilizing current Good Manufacturing Practices (cGMP) as overseen by the FDA and the International Standardization Organization (ISO) 13485:2016. These certifications attest to the quality standards set and achieved for the manufacture of our Drugs of Abuse products. Certified training is available on the Germaine Laboratories website for all registered customers.

THE TESTS YOU WANT. THE ACCURACY YOU NEED!



www.germainelabs.com

