



**GERMAINE**<sup>®</sup>  
LABORATORIES, INC.

# PRODUCT CATALOG 2022





## MISSION STATEMENT

**Quality In Everything We Do!** Through our actions, we earn the trust of our customers and team members. As a team, we are focused on meeting our customers' needs and keeping them for life. We believe continuous learning improves our contributions. Our work creates profit for our company and opportunity for us by creating quality and value for our customers.

## QUALITY PRODUCTS

All products sold by Germaine® Laboratories, Inc. are manufactured under strict quality control through cGMP (Good Manufacturing Practices). Our facility is in compliance with the Quality Management System Standards of ISO 13485:2016 and is applicable to the manufacture of In-Vitro Diagnostic Tests, In-Vitro Diagnostic Analyzers, and Reagents for Professional, Point of Care, and Self-Test Use used in the diagnosis, management, disease status and detection of Blood Glucose Monitoring, Cholesterol Monitoring, Clinical Chemistry, Drugs of Abuse, Hematology, Infectious Diseases, Pregnancy/Fertility, Serology, and Urinalysis.

## SPECIALIZED SERVICES

- ⇒ Bulk and OEM Packaging
- ⇒ Custom Labeling and Packaging
- ⇒ Custom Molded Cassettes
- ⇒ Regulatory consultation in FDA cGMP and ISO 13485:2016
- ⇒ Clinical Trials for Medical Devices and 510(k) Submissions
- ⇒ Contract Tableting Manufacturing



## TABLE OF CONTENTS

Flu A + B Tests..... 4

Blood Glucose Tests & Supplies ..... 5

Drugs Tests - Single Analytes ..... 6

Drugs Tests - MultiDrugs & Controls..... 7

Drug Cups & Collection Cups ..... 8

Febrile Antigens & Antisera ..... 9

Hemoglobin Tests & Supplies ..... 10

AimStrip® Tandem Measuring System ..... 11

Occult Blood FOB/GUAIAC ..... 12

Occult Blood IFOB ..... 13

Pregnancy Tests ..... 14

Ovulation Tests ..... 14

Thyroid Screen ..... 14

StrepAim® Strep A..... 15

AimStep® H.Pylori Tests ..... 15

AimStep® Mono Tests ..... 15

Serology Tests..... 15

Tablet Tests ..... 16

Tablet Contract Manufacturing ..... 16

Urinalysis Strips ..... 17

Microalbumin Strips ..... 17

AimStrip® Urine Analyzer 2..... 18

UA Controls & Supplies..... 18

VQPet H Blood Glucose Tests & Supplies ..... 19

COVID-19 Tests ..... 20



**CLIA  
WAIVED  
FOR SWAB  
SPECIMENS**

# FLU A+B TESTS



STORE AT 2° to 30°C

### SIMPLE PROCEDURE FOR SWAB SPECIMENS



#### BioSign® Flu A + B

**SIZE**      **REF#**  
25 tests      BSP-510-25

- ⇒ Rapid qualitative test that detects influenza type A & B
- ⇒ Includes Positive and Negative control swabs

### FLOCKED NASAL SWABS FOR IMPROVED PATIENT COMFORT



#### Flu A + B CONTROL SET PACKS

**REF#**

⇒ **2 Pack** .....(1 Negative • 1 Positive)      51002

⇒ **10 Pack** .....(5 Negative • 5 Positive)      51010




FOR  
SELF-TESTING  
AND  
PROFESSIONAL  
USE



STORE AT 2° to 30°C

# BLOOD GLUCOSE TESTS

Meter used for AimStrip® Plus Blood Glucose Strips | 5 YEAR WARRANTY



## AimStrip® Plus Blood Glucose Meter

REF#

37321

### ⇒ KIT INCLUDES:

Meter, 10 Lancets, Lancing Device, Carrying Case, User's Manual and Log Book

## AimStrip® Plus Blood Glucose Meter STARTER PACK

37355

### ⇒ KIT INCLUDES:

- 1 AimStrip® Plus Blood Glucose Meter (as above)
- Extras included: 1 AimStrip® Plus Test Strips (50/box) and 50 Lancets

Alternate Site Testing approved • English and Spanish instructions included • CLIA waived

CLIA  
WAIVED

## Using AimStrip® Plus Blood Glucose Meter



## AimStrip® Plus Blood Glucose Test Strips

SIZE REF#  
50 Strips 37350

- ⇒ 2 Bottles x 25 Strips per box
- ⇒ Code Chip included, does not change between lot numbers



## AimStrip® Plus Glucose Control Solution

SIZE REF#  
2 x 2mL 37312

- ⇒ Level 1 Control Solution
- ⇒ Level 2 Control Solution



## AimStrip® Plus Lancets

SIZE REF#  
100/box 37311

- ⇒ Sterile Lancets
- ⇒ 4 Bags x 25 each

AimStrip® Plus Lancing Device REF# 37331





# DRUGS TESTS

## SINGLE ANALYTES



STORE AT 2° to 30°C

PACKAGED IN BOXES OF 25 INDIVIDUAL TESTS

### ECONOMICAL URINE DIPDEVICE




	REF#		REF#
<b>AimScreen™ AMP</b> ⇒ for Amphetamine	41125	<b>AimScreen™ MTD</b> ⇒ for Methadone	46125
<b>AimScreen™ BAR</b> ⇒ for Barbiturate	42125	<b>AimScreen™ OPI</b> ⇒ for Opiates/Morphine, 2,000ng/mL	44125
<b>AimScreen™ BUP</b> ⇒ for Buprenorphine	42025	<b>AimScreen™ OXY</b> ⇒ for Oxycodone	49425
<b>AimScreen™ BZO</b> ⇒ for Benzodiazepine	42325	<b>AimScreen™ PCP</b> ⇒ for Phencyclidine	49125
<b>AimScreen™ COC</b> ⇒ for Cocaine, 300ng/mL	43125	<b>AimScreen™ PPX</b> ⇒ for Propoxyphene	49525
<b>AimScreen™ MET</b> ⇒ for Methamphetamine	47125	<b>AimScreen™ TCA</b> ⇒ for Tricyclic Antidepressants	44225
<b>AimScreen™ MDMA</b> ⇒ for Ecstasy	48225	<b>AimScreen™ THC</b> ⇒ for Marijuana	45125
<b>AimScreen™ MOP</b> ⇒ for Morphine/Opiates, 300ng/ML	48125		

### URINE CASSETTE | NOT CLIA WAIVED

	SIZE	REF#	SIZE	REF#
<b>AimStep® Cotinine</b> ⇒ for Cotinine	25 tests	40014	<b>AimScreen™ FEN</b> ⇒ for Fentanyl	42335



# DRUGS TESTS MULTIDRUGS



STORE AT 2° to 30°C

PACKAGED IN BOXES OF 25 INDIVIDUAL TESTS

## ECONOMICAL URINE DIPDEVICE



	REF#
<b>AimScreen™ MultiDrug 2</b> ⇒ for COC, THC	40225
<b>AimScreen™ MultiDrug 3</b> ⇒ for COC, MET, THC	40325
<b>AimScreen™ MultiDrug 4</b> ⇒ for COC, MET, OPI, THC	40425
<b>AimScreen™ MultiDrug 5</b> ⇒ for COC, AMP, MET, THC, OPI	45625
<b>AimScreen™ MultiDrug N5</b> ⇒ for COC, AMP, THC, OPI, PCP	40525
<b>AimScreen™ MultiDrug 6</b> ⇒ for COC, AMP, MET, THC, OPI, BZO	46325

	REF#
<b>AimScreen™ MultiDrug 6</b> ⇒ for COC, AMP, MET, THC, OPI, OXY	40625
<b>AimScreen™ MultiDrug 9</b> ⇒ for COC, AMP, MET, THC, OPI, PCP, BAR, BZO, MTD	40925
<b>AimScreen™ MultiDrug 10</b> ⇒ for COC, AMP, MET, THC, OPI, PCP, BZO, BAR, MTD, MDMA	40125
<b>AimScreen™ MultiDrug 12</b> ⇒ for COC, AMP, MET, THC, MTD, MDMA, OPI, OXY, PPX, PCP, BAR, BZO	40135

OTHER COMBINATIONS ARE AVAILABLE

## CONTROLS



STORE AT 2° to 8°C



	SIZE	REF#
<b>DOA Positive Control</b> <i>Positive Control (AMP, BAR, BUP, BZO, COC, MDMA, MET, MOP, MTD, OPI, OXY, PCP, PPX, TCA, THC)</i>	1 x 5mL	40002



<b>DOA Negative Control</b> <i>Negative control for all analytes</i>	1 x 5mL	40003
---	---------	-------



STORE AT 2° to 30°C


	SIZE	REF#
<b>AimStrip® Adulteration 7</b>	25 tests	40007

*Urine dipsticks for Creatinine, Nitrite, Glutaraldehyde, pH, Specific Gravity, Oxidant, PCC*






EXAMPLE OF VALID RESULTS  
5 PANEL CUP



EXAMPLE OF VALID RESULTS  
10 PANEL CUP




# SafeCup™ III DRUG CUPS



STORE AT 2° to 30°C

25 CUPS PER CASE

## URINE ROUND CUP



### SafeCup™ III 5 Panel

⇒ for COC, AMP, MET, THC, OPI

### SafeCup™ III N5 Panel

⇒ for COC, AMP, THC, OPI, PCP

### SafeCup™ III 6 Panel

⇒ for COC, AMP, THC, OPI, MET, BZO

### SafeCup™ III 10 Panel

⇒ for AMP, BAR, BZO, COC, MTD, MDMA, MET, OPI, PCP, THC

### SafeCup™ III 12 Panel

⇒ for AMP, BAR, BUP, BZO, COC, MTD, MDMA, MET, OPI, OXY, PCP, THC

### SafeCup™ III 12-M Panel

⇒ for COC, AMP, MET, THC, MOP, MTD, MDMA, OXY, PCP, BAR, BZO, BUP

### SafeCup™ III 12-MAD Panel

⇒ for COC, AMP, MET, THC, MOP, MTD, MDMA, OXY, PCP, BAR, BZO, BUP  
⇒ Includes 3 Adulterants OX/SG/PH

### SafeCup™ III 14 Panel

⇒ for COC, AMP, MET, THC, MOP, MTD, MDMA, OXY, PCP, BAR, BZO, BUP, TCA, PPX

REF#

45475

45485

40675

40105

41513

41577

41578

44418

OTHER COMBINATIONS ARE AVAILABLE

## COLLECTION CUPS



INSTRUCTIONS IN  
ENGLISH & SPANISH

### EasyAim™ Midstream Collection

SIZE

25

REF#

55525

- ⇒ Extra long handle allows easy collection
- ⇒ Includes: 1 Collection Cup, 1 extra long handle and 3 Towelettes

### Specimen Collection Cup

SIZE

400

REF#

182-030-1110

- ⇒ 90 x 53mm
- ⇒ Specimen Container with Green Cap NS with Temperature Strip



# FEBRILE ANTIGENS & ANTISERA



STORE AT 2° to 8°C

## BACTERIAL ANTIGENS

For testing patients serum using  
SLIDE OR TUBE TEST PROCEDURES



## CONTROLS

	SIZE	REF#
<b>Brucella Abortus Antiserum</b>	3mL	1100
⇒ for Brucella control		
<b>Febrile Positive Control</b>	5mL	1600
⇒ Positive control for B. abortus, P. Ox19, Sal. H a, Sal. H b, Sal. H d, Sal. O Gr. D		
<b>Febrile Negative Control</b>	5mL	1601
⇒ for all febrile antigens		
<b>Francisella Tularensis Antiserum</b>	2mL	1702
⇒ for F. tularensis control		
<b>Proteus Poly Antiserum</b>	3mL	1303
⇒ Positive control for Proteus OX2, OX19, OXK		
<b>Salmonella H Poly Positive Control</b>	3mL	1606
⇒ Positive control for Sal. H a,b,c,d		
<b>Salmonella O Poly Positive Control</b>	3mL	1611
⇒ Positive control Sal. O Gr. A, B, C, D		

	SIZE	REF#
<b>Brucella Abortus Antigen</b>	5mL	1000
<b>Brucella Melitensis Antigen</b>	5mL	1002
<b>Francisella Tularensis Antigen</b>	3mL	1703
<b>Proteus OX2 Antigen</b>	5mL	1201
<b>Proteus OX19 Antigen</b>	5mL	1202
<b>Proteus OXK Antigen</b>	5mL	1200
<b>Salmonella H a Antigen (paratyphoid A)</b>	5mL	1400
<b>Salmonella H b Antigen (paratyphoid B)</b>	5mL	1401
<b>Salmonella H c Antigen (paratyphoid C)</b>	5mL	1402
<b>Salmonella H d Antigen (typhoid H)</b>	5mL	1403
<b>Salmonella O Gr. A Antigen (somatic ag.)</b>	5mL	1500
<b>Salmonella O Gr. B Antigen (somatic ag.)</b>	5mL	1501
<b>Salmonella O Gr. C Antigen (somatic ag.)</b>	5mL	1502
<b>Salmonella O Gr. D Antigen (typhoid O)</b>	5mL	1503

CLIA  
WAIVED

# HEMOGLOBIN TESTS & SUPPLIES



STORE AT 2° to 30°C

STARTER  
KIT 200



## RAPID HEMOGLOBIN & HEMATOCRIT RESULTS

### AimStrip® Hemoglobin Starter Kit 200

SIZE  
200 tests REF#  
78200

- ⇒ 1 AimStrip® Hemoglobin Meter
- ⇒ 200 Test Cartridges
- ⇒ 2 Optical Verifiers
- ⇒ 1 Carrying Case
- ⇒ User's Manual
- ⇒ 240 Capillary Tubes
- ⇒ 4 Code Chips
- ⇒ 4 AAA Batteries
- ⇒ 4 Results Recording Pads



STORE AT 2° to 30°C



### AimStrip® Hemoglobin Meter

SIZE REF#  
1 meter 78101

- ⇒ Includes:  
1 Meter, 2 Optical Verifiers,  
1 Carrying Case, 4 AAA Batteries



### AimStrip® Hemoglobin Test Combo

SIZE REF#  
50/60 combo 78250

- ⇒ 50 Test Cartridges
- ⇒ 60 Capillary Tubes
- ⇒ CLIA Waived



STORE AT 2° to 8°C



### AimStrip® Hb Control Set

SIZE REF#  
3 x 2mL 78103

- ⇒ 1 x 2mL Control Solution 0 (low)
- ⇒ 1 x 2mL Control Solution 1 (normal)
- ⇒ 1 x 2mL Control Solution 2 (high)



### AimStrip® Hb Optical Verifiers

SIZE REF#  
2 pack 78102

- ⇒ 2 Optical Verifiers



### AimStrip® Capillary Transfer Tubes

SIZE REF#  
50 pack 78155

- ⇒ 50 Capillary Tubes



### AimStrip® Hb AC Adaptor

SIZE REF#  
1 each 78111

- ⇒ 1 AC Adaptor (for meter)

CLIA  
WAIVED



REF# 77301

# AimStrip® TANDEM MEASURING SYSTEM

## LIPID PROFILE & GLUCOSE MEASURING SYSTEM

### METER USED FOR TESTING:

- AimStrip® Total Cholesterol Test Strips
- AimStrip® Triglycerides Test Strips
- AimStrip® Blood Glucose Test Strips

AimStrip® TANDEM BRAND  
TEST STRIPS & CONTROL  
SOLUTIONS ARE REQUIRED  
BUT NOT INCLUDED IN KIT.



STORE AT 2° to 30°C



MULTIPLE PATIENT USE



1 KIT

INCLUDES:

- ➔ AimStrip® Tandem Meter
- ➔ Operation Manual
- ➔ Warranty Card
- ➔ Carrying Case
- ➔ 2 AAA Batteries
- ➔ Log Book



## AimStrip® Tandem TEST STRIPS - For use with the AimStrip® Tandem Measuring System



**AimStrip® Tandem  
TOTAL CHOLESTEROL TEST STRIPS**

SIZE	REF#
10 Strips	77310

- ➔ 10 Test Strips
- ➔ 10 Capillary Rods



**AimStrip® Tandem  
TRIGLYCERIDES TEST STRIPS**

SIZE	REF#
10 Strips	77410

- ➔ 10 Test Strips
- ➔ 10 Capillary Rods



**AimStrip® Tandem  
BLOOD GLUCOSE TEST STRIPS**

SIZE	REF#
50 Strips	77550

- ➔ Packed in 2 Bottles
- ➔ 25 Test Strips per Bottle

## AimStrip® Tandem CONTROL SOLUTIONS

CONTROLS ARE FOR USE WITH AimStrip® TANDEM METER ONLY.

**AimStrip® Tandem TOTAL  
CHOLESTEROL Control Solution**

SIZE	REF#
2 x 3mL	77312

- ➔ Control Solution Level 1 (3mL)
- ➔ Control Solution Level 2 (3mL)

**AimStrip® Tandem  
TRIGLYCERIDES Control Solution**

SIZE	REF#
2 x 3mL	77412

- ➔ Control Solution Level 1 (3mL)
- ➔ Control Solution Level 2 (3mL)

**AimStrip® Tandem  
GLUCOSE Control Solution**

SIZE	REF#
3 x 3.5mL	77512

- ➔ Control Solution Level 1 (3.5mL)
- ➔ Control Solution Level 2 (3.5mL)
- ➔ Control Solution Level 3 (3.5mL)



STORE AT 8° to 30°C

CLIA  
WAIVED

# OCCULT BLOOD FOB/GUAIAC

**GUAIAC METHOD for detecting occult blood in stool**



STORE AT 2° to 30°C



**COMPLIANCE MONITORING**  
HIPAA Compliant program that reminds patients to return Collection Kits.

### COMPLIANCE VERSION:

- ⇒ Full compliance version includes Compliance Monitoring
- ⇒ Return Labels and "Postage Pre-Paid" on Return Mailers
- ⇒ HIPAA compliant Patient Agreement forms

### EACH KIT PROVIDES THE FOLLOWING:

- ⇒ On-Slide controls, detects Hb at 0.6mg of stool
- ⇒ Enhanced Readability results in 30 seconds
- ⇒ Collection Sticks & Developing Solution



### FULL COMPLIANCE VERSION

The following kit **INCLUDES Compliance Monitoring and Pre-Paid Postage.**

**Compliance Gold® FOB Test-ER, TRI-SLIDES**

REF#  
74252

### ECONOMY VERSION

The following kits **DO NOT include Compliance Monitoring and Postage.**

**Compliance Gold® FOB Test-ER, TRI-SLIDES**

REF#  
74250

**Compliance Gold® FOB Test-ER, SINGLES**

REF#  
74200

### DEVELOPER

	SIZE	REF#
<b>Compliance Gold® FOB Developer</b>	8 x 10 mL	74208

**TRI-SLIDES:  
50 TESTS PER KIT**



**SINGLES:  
100 TESTS PER KIT**



ENGLISH  
AND  
SPANISH  
INSTRUCTIONS  
AVAILABLE



**CLIA  
WAIVED**

# OCCULT BLOOD iFOB



STORE AT 2° to 30°C

30 TESTS PER KIT



**IMMUNOLOGICAL TEST for occult blood in stool specimens**



**COMPLIANCE MONITORING**  
HIPAA Compliant program that reminds patients to return Collection Kits.

**COMPLIANCE VERSION:**

- ⇒ Full compliance version includes Compliance Monitoring
- ⇒ Return Labels & "Postage Pre-Paid" on Return Mailers
- ⇒ HIPAA compliant Patient Agreement forms



**FULL COMPLIANCE VERSION**

The following kit **INCLUDES Compliance Monitoring and Pre-Paid Postage.**

**Compliance Gold® iFOB** REF# 74034

- ⇒ Includes Test Devices, Collection Bottles & Collection Papers

**ECONOMY VERSION**

The following kit **DOES NOT include Compliance Monitoring and Postage.** (Return Mailers without postage are included).

**Compliance Gold® iFOB Simple** REF# 74030

- ⇒ Includes Test Devices, Collection Bottles & Collection Papers



STORE AT 2° to 30°C



- 15 Sets:
- ⇒ Return Mailers
  - ⇒ Collection Bottles
  - ⇒ Collection Papers

**Patient Collection Kits** REF# 74014  
**With Pre-Paid Postage**

*Use with Compliance Gold® iFOB*  
**Without Postage** 74015

*Use with Compliance Gold® iFOB Simple*



**Compliance Gold® iFOB Easy** REF# 74033  
**SIZE** 30 tests

- ⇒ Test Devices
- ⇒ Collection Bottles
- ⇒ No Compliance, Mailers or Postage



STORE AT 2° to 8°C



**iFOB Controls** REF# 74211  
**SIZE** 2 x 1mL

- ⇒ Controls for occult blood tests
- ⇒ 1 x 1mL (Positive Control)
- ⇒ 1 x 1mL (Negative Control)

# PREGNANCY TESTS



**CLIA  
WAIVED  
FOR  
URINE**



STORE AT 2° to 30°C



**AimStep® Pregnancy**  
**SIZE**      **REF#**  
 30 Tests    97030

- ⇒ Urine cassette method
- ⇒ Sensitive to 20 mIU/mL hCG
- ⇒ CLIA waived, also available in OTC



**AimStep® Combo Pregnancy**  
**SIZE**      **REF#**  
 30 Tests    97730

- ⇒ Serum / Urine cassette method
- ⇒ Sensitive to 20 mIU/mL hCG
- ⇒ CLIA waived for urine testing



**AimStrip® Pregnancy**  
**SIZE**      **REF#**  
 30 Tests    88530

- ⇒ Urine dipstick method
- ⇒ Sensitive to 20 mIU/mL hCG
- ⇒ CLIA waived, also available in OTC



STORE AT 2° to 30°C

## AccuHome® PREGNANCY

**OTC URINE MIDSTREAM METHOD**



**SIZE**      **REF#**  
 1 Test      89101  
 2 Tests    89102

- ⇒ OTC final box packing
- ⇒ Sensitive to 25 mIU/mL hCG
- ⇒ CLIA waived



STORE AT 2° to 8°C

## HCG CONTROLS

FOR USE WITH URINE PREGNANCY TESTS



**HCG Urine Controls**  
**SIZE**      **REF#**  
 2 x 5mL    71002

- ⇒ Negative Control (1 x 5mL)
- ⇒ Positive Control (1 x 5mL)

FOR USE WITH SERUM PREGNANCY TESTS



**HCG Serum Controls**  
**SIZE**      **REF#**  
 3 x 2.5mL    71003

- ⇒ Negative, Low Positive, and High Positive Control

**CLIA  
WAIVED**

## OVULATION TESTS



**AimStep® Ovulation**  
**SIZE**      **REF#**  
 30 Tests    75030

- ⇒ Urine cassette method
- ⇒ OTC or professional use



**AccuHome® Ovulation**  
**SIZE**      **REF#**  
 1 Test      75501

- ⇒ OTC Urine midstream method
- ⇒ Bulk packed

**CLIA  
WAIVED**

## THYROID SCREEN



**AimStep® Thyroid Screen**  
**SIZE**      **REF#**  
 20 Tests    76220

- ⇒ Fingerstick Whole Blood cassette method
- ⇒ Detects TSH, Professional use

DIPSTICK METHOD FOR RAPID  
DETECTION OF GROUP A STREP

# StrepAim® STREP A



**CLIA  
WAIVED**



- KITS INCLUDE:
- Dipsticks
  - Controls
  - Extract Reagents
  - Swabs
  - Tubes
  - Workstation

TEST DIRECTLY FROM  
THROAT SWABS



STORE AT 2° to 30°C



	SIZE	REF#
StrepAim®	40 Tests	73040
StrepAim®	25 Tests	73025



STORE AT 2° to 30°C



SERUM, PLASMA AND WHOLE BLOOD CASSETTE METHOD

# AimStep® H. PYLORI

# AimStep® MONO

**CLIA  
WAIVED  
FOR WHOLE  
BLOOD**



SIZE REF#  
20 tests 33070



SIZE REF#  
25 tests 77025

KITS INCLUDE:



➤ H.Pylori Cassettes, Positive and Negative Controls, Developer Solution & Fingersticks

KITS INCLUDE:



➤ Mono Cassettes, Positive and Negative Controls, Developer Solution & Fingersticks



STORE AT 2° to 8°C

# SEROLOGY



QUALITATIVE & SEMI-QUANTITATIVE PROCEDURES

**Latex agglutination SLIDE METHOD**

KITS INCLUDE: Latex, Controls, Buffer, Test Slide, and Pipettes

SIZE REF#  
**Aim™ ASO Test** 50 tests 31050

➤ for non-diluted serum to detect for anti-streptolysin O

SIZE REF#  
**Aim™ CRP Test** 100 tests 32100

➤ for non-diluted serum to detect for C-Reactive Protein

SIZE REF#  
**Aim™ RF Test** 100 tests 35100

➤ for non-diluted serum to detect for Rheumatoid Factor

**Carbon flocculation SLIDE METHOD**

KIT INCLUDES: Carbon Antigen, Controls, Disposable Slides, and Pipettes

SIZE REF#  
**Aim™ RPR Test** 500 tests 36500

➤ for serum or plasma to detect for reagin antibodies

# TABLET TESTS



STORE AT 2° to 30°C

## AimTab™ BILIRUBIN TABLETS

REF# 13200  
SIZE 100 tests

Tablet method for Bilirubin in urine.

- ⇒ Compare to Ictotest®, trademark by Siemens Healthcare
- ⇒ Test Mats included
- ⇒ Controls, see UA Liquid Controls Plus, Ref #51116



STORE AT 10° to 35°C

## AimTab™ BLUING TABLETS

REF# 40021  
SIZE 100 tests



Tablet method designed for use in all drug testing collections to prevent adulteration.

- ⇒ Meets Federal DOT Regulation 49 CFR 40.41
- ⇒ Not classified as an In-vitro Diagnostic test



STORE AT 2° to 30°C



## AimTab™ KETONE TABLETS

Tablet method for the semi-quantitative determination of ketones (acetoacetic acid and acetone) in urine, serum, plasma & whole blood.

- ⇒ Detects as little as 5 mg/dL of acetoacetic acid in urine
- ⇒ Compare to Acetest®, trademark by Bayer Healthcare LLC.
- ⇒ Urine Control, see UA Liquid Controls Plus, Ref #51116
- ⇒ CLIA moderately complex (testing serum, plasma or whole blood)



**Ketone Serum Controls**  
SIZE 2 x 2mL  
REF# 13112

- ⇒ 1 Negative Control
- ⇒ 1 Positive Control
- ⇒ **Store at 2° to 8°C**

## AimTab™ REDUCING SUBSTANCES TABLETS



STORE AT 2° to 30°C

Tablet method for the semi-quantitative determination of the amount of reducing substances (glucose, lactose, fructose, galactose, pentoses) in urine, stool, or wine.

- ⇒ Wine making, used to estimate the residual sugar level of wine
- ⇒ Compare to Clinitest®, trademark by Bayer Healthcare LLC.
- ⇒ For Forensic Use Only



	SIZE	REF#
<b>BOTTLED TABLETS</b>	36 tests	13336
<b>POUCHED TABLETS</b>	20 tests	13320


- ⇒ 1 tablet per sealed foil pouch

**TABLETING CONTRACT MANUFACTURING**

A small-scale, personal contract tableting manufacturer. Focused on full development support for customers with small-scale manufacturing needs. cGMP compliant contract manufacturing services are available for any formulation. Compliant with Title 21 CFR part 11. CONTACT CUSTOMER SERVICE.



# AimStrip® URINALYSIS STRIPS



STORE AT 2° to 30°C

PACKAGED IN BOTTLES  
OF 100 REAGENT STRIPS

## URINE DIPSTICKS

	REF#		REF#
<b>AimStrip® US-G</b>	50100	<b>AimStrip® US-5</b>	55100
⇒ for Glucose		⇒ for Glucose, Protein, Ketone, pH, Blood	
<b>AimStrip® US-K</b>	50121	<b>AimStrip® US-5B</b>	55200
⇒ for Ketone		⇒ for Blood, Glucose, Protein, Nitrite, Leukocytes	
<b>AimStrip® US-2</b>	52100	<b>AimStrip® US-9</b>	59100
⇒ for Glucose & Protein		⇒ for Glucose, Bilirubin, Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite	
<b>AimStrip® US-3</b>	53100	<b>AimStrip® 10-SG</b>	51100
⇒ for Glucose, Protein, pH		⇒ for Glucose, Bilirubin, Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite, Leukocytes	
<b>AimStrip® US-4OB</b>	54100		
⇒ for Glucose, Protein, Nitrite, Leukocytes			

PLEASE INQUIRE ABOUT OTHER AVAILABLE COMBINATIONS

## UTI COMBO STRIPS



SIZE REF#  
1 Test 52301

- ⇒ for Leukocytes and Nitrites
- ⇒ For OTC use

## MICROALBUMIN STRIPS



### AimStrip® MACRE Strips

SIZE REF#  
25 Tests 52275

- ⇒ Sensitivity as low as 2.0 mg/dL
- ⇒ Detects Microalbumin & Creatinine as low as 10mg/L
- ⇒ Semi-Quantitative A:C ratio results in one minute
- ⇒ Same procedure as leading national brand



EXCLUSIVELY  
READ ON  
AimStrip® URINE  
ANALYZER 2



### MicroAlbumin 2-1 Combo Strips

25 Tests 52255

- ⇒ Urine Dipstick for Microalbumin and Creatinine
- ⇒ Determines Albumin to Creatinine ratio
- ⇒ For visual read only, analyzer not required

FOR THE SEMI-QUANTITATIVE DETECTION  
OF ANALYTES IN URINE



# AimStrip® URINE ANALYZER 2

REF# 51222



STORE AT 2° to 30°C



1 UNIT

- ⇒ Reads all combinations of strips
- ⇒ Single Test Mode, reads 1 strip at a time
- ⇒ Continuous Test Mode, reads 120 strips/hour
- ⇒ Automatic calibration for accurate results
- ⇒ Automatically saves the last 2000 results
- ⇒ Includes RS232C and USB ports to export results

LANGUAGES INCLUDED: ENGLISH, FRENCH, PORTUGUESE, SPANISH, ITALIAN



REF# 51110  
Barcode reader included

# AimStrip® URINE AUTO ANALYZER

LANGUAGES INCLUDED: ENGLISH, SPANISH, PORTUGUESE

- ⇒ High throughput, up to 500 tests per hour
- ⇒ Large touch screen for easy navigation
- ⇒ Reads 4, 8, 10 parameter strips
- ⇒ Automatic calibration for accurate results
- ⇒ Automatically saves the last 2000 results
- ⇒ Includes RS232C and USB ports to export results



STORE AT 2° to 8°C

## UA CONTROLS

### Liquid



	SIZE	REF#
<b>UA Liquid Controls</b>	2 x 5mL	51115
⇒ Negative Control (1 x 5mL), Positive Control (1 x 5mL)		
<b>UA Liquid Controls</b>	6 x 10mL	51161
⇒ Negative Control (3 x 10mL), Positive Control (3 x 10mL)		
<b>UA Liquid Controls Plus</b>	2 x 5mL	51116
⇒ Negative Control (1 x 5mL), Positive Control (1 x 5mL)		

### DipTube



	SIZE	REF#
<b>UA DipTube Controls</b>	2 x 12mL	51112
⇒ Negative Control (1 x 12mL), Positive Control (1 x 12mL)		
<b>UA DipTube Controls</b>	4 x 12 mL	51142
⇒ Negative Control (2 x 12mL), Positive Control (2 x 12mL)		



STORE AT 2° to 30°C

## PAPER

4 ROLLS/PK



	REF#
<b>Printer Paper</b>	51224
⇒ Thermal paper to print 200 results per roll	
<b>Printer Peel &amp; Stick Paper</b>	51234
⇒ Peel & stick printer paper to print 100 results/roll	
⇒ Conveniently place results into patient records	

## BARCODE READER



	REF#
<b>Barcode Reader For Analyzers</b>	51320
⇒ Use with AimStrip® Urine Analyzer 2 & AimStrip® Urine Auto Analyzer to record patient ID	



FOR  
PROFESSIONAL  
USE  
ONLY  
FOR PETS



# VQPET H BLOOD GLUCOSE TESTS & SUPPLIES

FRESH WHOLE BLOOD (CAPILLARY, VENOUS, ARTERIAL)  
CAN BE USED FOR TESTING FOR VETERINARY CLINICS.



CALIBRATED FOR  
DOGS & CATS



GLUCOSE  
DEHYDROGENASE-  
FLAVIN ADENINE  
DINUCLEOTIDE  
(GDH-FAD) ENZYME  
USED



## VQPET H GLUCOSE STARTER KIT

REF# VGM30-431



### 1 KIT

#### INCLUDES:

- ⇒ 1 Pet Blood Glucose Meter
- ⇒ 10 Test Strips
- ⇒ 10 Sterile Lancets
- ⇒ 1 Lancing Device
- ⇒ 2mL Control Solution
- ⇒ Carrying Case
- ⇒ User's Manual
- ⇒ Quick Reference Guide
- ⇒ Logbook

## TEST STRIPS



VQPET  
GLUCOSE TEST STRIPS

REF#  
VGS01-421

- ⇒ 50 Strips/Box
- ⇒ (2) 25 Strips/Vial

## CONTROL SOLUTION

VQPET  
CONTROL SOLUTION



REF#  
VGM30-371

- ⇒ 1 x 2mL

Support Alternative Site Testing for Pets • Only 0.5µL Blood Sample Required from Pet



# COVID-19 TESTS



STORE AT 2° to 30°C

QUALITATIVE & RAPID DETECTION OF CORONAVIRUS INFECTION

## Status™ COVID-19/FLU A&B (SWAB)

REF# 33225  
SIZE 25 tests

- ⇒ 1 test + 1 specimen = 3 Results
- ⇒ Rapid Immunoassay
- ⇒ Flocked Nasopharyngeal Swab
- ⇒ COVID-19 Sensitivity: 93.9%
- ⇒ COVID-19 Specificity: 100%
- ⇒ FDA Emergency Use Authorized (EUA)



- KITS INCLUDE:
- ⇒ 25 Test Devices
  - ⇒ 25 Extraction Reagent in Capsules
  - ⇒ 25 Sterile Swabs
  - ⇒ 1 Positive Control Swab
  - ⇒ 1 Negative Control Swab
  - ⇒ 1 Instructions for Use
  - ⇒ 1 Quick Reference Instructions

## Celltrion DiaTrust™ COVID-19 AG RAPID TEST (SWAB)

REF# CT-P60 D-2 02  
SIZE 25 tests

- ⇒ Lateral Flow Immunoassay
- ⇒ Nasopharyngeal Swab (NPS)
- ⇒ Clinical Sensitivity: 93.33% (28/30)
- ⇒ Clinical Specificity: 99.03% (102/103)
- ⇒ FDA Emergency Use Authorized (EUA)



- KITS INCLUDE:
- ⇒ 25 Test Devices
  - ⇒ 25 Sterile Swabs
  - ⇒ 25 Extraction Buffer
  - ⇒ 25 Filter Caps
  - ⇒ 1 Positive Control Swab
  - ⇒ 1 Negative Control Swab
  - ⇒ 1 Quick Reference Instructions



## FaStep ANTIBODY RAPID DIAGNOSTIC TEST

REF# COV-W23  
SIZE 20 tests

- ⇒ Combined Sensitivity: 100%
- ⇒ Combined Specificity: 98.8%
- ⇒ FDA Emergency Use Authorized (EUA)



- KITS INCLUDE:
- ⇒ 20 Tests
  - ⇒ 20 Disposable Pipettes
  - ⇒ 20 Sterile Safety Lancets
  - ⇒ 20 Alcohol Prep Pads
  - ⇒ Buffer
  - ⇒ Package Insert





[www.germainelabs.com](http://www.germainelabs.com)

## FOR FURTHER INFORMATION CONTACT CUSTOMER SERVICE



800.854.8446  
210.692.4192  
210.692.4198 Fax

**Hours** ..... 8:30am to 5:00pm  
Central Standard Time

**Errors** ..... Products received in error  
must be reported within 10 days

**Damaged** ..... Products damaged in transit  
must be filed within 10 days

**P.O.D** ..... Requests for proof of delivery  
must be made within 60 days

**Returns** ..... All returns must be authorized.  
Products not returnable include  
expired or short dated products,  
opened products, and products  
improperly handled or stored.

**Remit to:** ..... PO Box 47015, Dept. R  
San Antonio, Texas 78265

**Shipping Address:** ..... 11030 Wye Drive  
San Antonio, Texas 78217 USA



**THE TESTS YOU WANT. THE ACCURACY YOU NEED!**