



PRODUCT CATALOG 2023







MISSION STATEMENT

Quality In Everything We Do! Through our actions, we earn the trust of our customers and team members. As a team, we are focused on meeting our customers' needs and keeping them for life. We believe continuous learning improves our contributions. Our work creates profit for our company and opportunity for us by creating quality and value for our customers.

QUALITY PRODUCTS

All products sold by Germaine® Laboratories, Inc. are manufactured under strict quality control through cGMP (Good Manufacturing Practices). Our facility is in compliance with the Quality Management System Standards of ISO 13485:2016 and is applicable to the manufacture of In-Vitro Diagnostic Tests, In-Vitro Diagnostic Analyzers, and Reagents for Professional, Point of Care, and Self-Test Use used in the diagnosis, management, disease status and detection of Blood Glucose Monitoring, Cholesterol Monitoring, Clinical Chemistry, Drugs of Abuse, Hematology, Infectious Diseases, Pregnancy/Fertility, Serology, and Urinalysis.

SPECIALIZED SERVICES

- Bulk and OEM Packaging
- Custom Labeling and Packaging
- Custom Molded Cassettes
- Regulatory consultation in FDA cGMP and ISO 13485:2016
- Clinical Trials for Medical Devices and 510(k) Submissions
- Contract Tableting Manufacturing



TABLE OF CONTENTS

Flu A + B Tests..... 4

Blood Glucose Tests & Supplies 5

Drugs Tests - Single Analytes 6

Drugs Tests - MultiDrugs & Controls..... 7

Drug Cups & Collection Cups 8

Febrile Antigens & Antisera 9

Hemoglobin Tests & Supplies 10

AimStrip® Tandem Measuring System 11

Occult Blood FOB/GUAIAC 12

Occult Blood IFOB 13

Pregnancy Tests 14

Ovulation Tests 14

Thyroid Screen 14

StrepAim® Strep A..... 15

AimStep® H.Pylori Tests 15

AimStep® Mono Tests 15

Serology Tests..... 15

Tablet Tests 16

Tablet Contract Manufacturing 16

Urinalysis Strips 17

Microalbumin Strips 17

AimStrip® Urine Analyzer 2..... 18

UA Controls & Supplies..... 18

VQPet H Blood Glucose Tests & Supplies 19

COVID-19 Tests 20



**CLIA
WAIVED
FOR SWAB
SPECIMENS**

FLU A+B TESTS



STORE AT 2° to 30°C

SIMPLE PROCEDURE FOR SWAB SPECIMENS



BioSign® Flu A + B

SIZE **REF#**
25 tests BSP-510-25

- ⇒ Rapid qualitative test that detects influenza type A & B
- ⇒ Includes Positive and Negative control swabs

FLOCKED NASAL SWABS FOR IMPROVED PATIENT COMFORT



Flu A + B CONTROL SET PACKS

REF#

⇒ **2 Pack**(1 Negative • 1 Positive) 51002

⇒ **10 Pack**(5 Negative • 5 Positive) 51010



FOR
SELF-TESTING
AND
PROFESSIONAL
USE



STORE AT 2° to 30°C

BLOOD GLUCOSE TESTS

Meter used for AimStrip® Plus Blood Glucose Strips | 5 YEAR WARRANTY



AimStrip® Plus Blood Glucose Meter

REF#

37321

⇒ KIT INCLUDES:

Meter, 10 Lancets, Lancing Device, Carrying Case, User's Manual and Log Book

AimStrip® Plus Blood Glucose Meter STARTER PACK

37355

⇒ KIT INCLUDES:

- 1 AimStrip® Plus Blood Glucose Meter (as above)
- Extras included: 1 AimStrip® Plus Test Strips (50/box) and 50 Lancets

Alternate Site Testing approved • English and Spanish instructions included • CLIA waived



Using AimStrip® Plus Blood Glucose Meter



AimStrip® Plus Blood Glucose Test Strips

SIZE REF#
50 Strips 37350

- ⇒ 2 Bottles x 25 Strips per box
- ⇒ Code Chip included, does not change between lot numbers



AimStrip® Plus Glucose Control Solution

SIZE REF#
2 x 2mL 37312

- ⇒ Level 1 Control Solution
- ⇒ Level 2 Control Solution



AimStrip® Plus Lancets

SIZE REF#
100/box 37311

- ⇒ Sterile Lancets
- ⇒ 4 Bags x 25 each

AimStrip® Plus Lancing Device REF# 37331





DRUGS TESTS

SINGLE ANALYTES



STORE AT 2° to 30°C

PACKAGED IN BOXES OF 25 INDIVIDUAL TESTS

ECONOMICAL URINE DIPDEVICE



	REF#		REF#
AimScreen™ AMP ⇒ for Amphetamine	41125	AimScreen™ MTD ⇒ for Methadone	46125
AimScreen™ BAR ⇒ for Barbiturate	42125	AimScreen™ OPI ⇒ for Opiates/Morphine, 2,000ng/mL	44125
AimScreen™ BUP ⇒ for Buprenorphine	42025	AimScreen™ OXY ⇒ for Oxycodone	49425
AimScreen™ BZO ⇒ for Benzodiazepine	42325	AimScreen™ PCP ⇒ for Phencyclidine	49125
AimScreen™ COC ⇒ for Cocaine, 300ng/mL	43125	AimScreen™ PPX ⇒ for Propoxyphene	49525
AimScreen™ MET ⇒ for Methamphetamine	47125	AimScreen™ TCA ⇒ for Tricyclic Antidepressants	44225
AimScreen™ MDMA ⇒ for Ecstasy	48225	AimScreen™ THC ⇒ for Marijuana	45125
AimScreen™ MOP ⇒ for Morphine/Opiates, 300ng/ML	48125		

URINE CASSETTE | NOT CLIA WAIVED


	SIZE	REF#
AimStep® Cotinine ⇒ for Cotinine	25 tests	40014
AimScreen™ FEN ⇒ for Fentanyl	25 tests	42335

STRIPS | NOT CLIA WAIVED

	SIZE	REF#
AimStrip FEN ⇒ for Fentanyl	25 tests	43345



DRUGS TESTS MULTIDRUGS



STORE AT 2° to 30°C

PACKAGED IN BOXES OF 25 INDIVIDUAL TESTS

ECONOMICAL URINE DIPDEVICE



	REF#
AimScreen™ MultiDrug 2 ⇒ for COC, THC	40225
AimScreen™ MultiDrug 3 ⇒ for COC, MET, THC	40325
AimScreen™ MultiDrug 4 ⇒ for COC, MET, OPI, THC	40425
AimScreen™ MultiDrug 5 ⇒ for COC, AMP, MET, THC, OPI	45625
AimScreen™ MultiDrug N5 ⇒ for COC, AMP, THC, OPI, PCP	40525
AimScreen™ MultiDrug 6 ⇒ for COC, AMP, MET, THC, OPI, BZO	46325

	REF#
AimScreen™ MultiDrug 6 ⇒ for COC, AMP, MET, THC, OPI, OXY	40625
AimScreen™ MultiDrug 9 ⇒ for COC, AMP, MET, THC, OPI, PCP, BAR, BZO, MTD	40925
AimScreen™ MultiDrug 10 ⇒ for COC, AMP, MET, THC, OPI, PCP, BZO, BAR, MTD, MDMA	40125
AimScreen™ MultiDrug 12 ⇒ for COC, AMP, MET, THC, MTD, MDMA, OPI, OXY, PPX, PCP, BAR, BZO	40135

OTHER COMBINATIONS ARE AVAILABLE

CONTROLS



STORE AT 2° to 8°C



	SIZE	REF#
DOA Positive Control <i>Positive Control (AMP, BAR, BUP, BZO, COC, MDMA, MET, MOP, MTD, OPI, OXY, PCP, PPX, TCA, THC)</i>	1 x 5mL	40002



DOA Negative Control <i>Negative control for all analytes</i>	1 x 5mL	40003
---	---------	-------



STORE AT 2° to 30°C


	SIZE	REF#
AimStrip® Adulteration 7	25 tests	40007

Urine dipsticks for Creatinine, Nitrite, Glutaraldehyde, pH, Specific Gravity, Oxidant, PCC






EXAMPLE OF VALID RESULTS
5 PANEL CUP



EXAMPLE OF VALID RESULTS
10 PANEL CUP



SafeCup™ III DRUG CUPS



STORE AT 2° to 30°C

25 CUPS PER CASE

URINE ROUND CUP



SafeCup™ III 5 Panel

⇒ for COC, AMP, MET, THC, OPI

SafeCup™ III N5 Panel

⇒ for COC, AMP, THC, OPI, PCP

SafeCup™ III 6 Panel

⇒ for COC, AMP, THC, OPI, MET, BZO

SafeCup™ III 10 Panel

⇒ for AMP, BAR, BZO, COC, MTD, MDMA, MET, OPI, PCP, THC

SafeCup™ III 12-M Panel

⇒ for COC, AMP, MET, THC, MOP, MTD, MDMA, OXY, PCP, BAR, BZO, BUP

SafeCup™ III 12-MAD Panel

⇒ for COC, AMP, MET, THC, MOP, MTD, MDMA, OXY, PCP, BAR, BZO, BUP

⇒ Includes 3 Adulterants OX/SG/PH

SafeCup™ III 14 Panel

⇒ for COC, AMP, MET, THC, MOP, MTD, MDMA, OXY, PCP, BAR, BZO, BUP, TCA, PPX

REF#

45475

45485

40675

40105

41577

41578

44418

OTHER COMBINATIONS ARE AVAILABLE

COLLECTION CUPS



INSTRUCTIONS IN
ENGLISH & SPANISH

EasyAim™ Midstream Collection

SIZE

25

REF#

55525

- ⇒ Extra long handle allows easy collection
- ⇒ Includes: 1 Collection Cup, 1 extra long handle and 3 Towelettes

Specimen Collection Cup

SIZE

400

REF#

182-030-1110

- ⇒ 90 x 53mm
- ⇒ Specimen Container with Green Cap NS with Temperature Strip



FEBRILE ANTIGENS & ANTISERA



STORE AT 2° to 8°C

BACTERIAL ANTIGENS

For testing patients serum using
SLIDE OR TUBE TEST PROCEDURES



CONTROLS

	SIZE	REF#
Brucella Abortus Antiserum	3mL	1100
⇒ for Brucella control		
Febrile Positive Control	5mL	1600
⇒ Positive control for B. abortus, P. Ox19, Sal. H a, Sal. H b, Sal. H d, Sal. O Gr. D		
Febrile Negative Control	5mL	1601
⇒ for all febrile antigens		
Francisella Tularensis Antiserum	2mL	1702
⇒ for F. tularensis control		
Proteus Poly Antiserum	3mL	1303
⇒ Positive control for Proteus OX2, OX19, OXK		

	SIZE	REF#
Brucella Abortus Antigen	5mL	1000
Brucella Melitensis Antigen	5mL	1002
Francisella Tularensis Antigen	3mL	1703
Proteus OX2 Antigen	5mL	1201
Proteus OX19 Antigen	5mL	1202
Proteus OXK Antigen	5mL	1200
Salmonella H a Antigen (paratyphoid A)	5mL	1400
Salmonella H b Antigen (paratyphoid B)	5mL	1401
Salmonella H c Antigen (paratyphoid C)	5mL	1402
Salmonella H d Antigen (typhoid H)	5mL	1403
Salmonella O Gr. A Antigen (somatic ag.)	5mL	1500
Salmonella O Gr. B Antigen (somatic ag.)	5mL	1501
Salmonella O Gr. C Antigen (somatic ag.)	5mL	1502
Salmonella O Gr. D Antigen (typhoid O)	5mL	1503



HEMOGLOBIN TESTS & SUPPLIES



STORE AT 2° to 30°C

STARTER KIT 200



RAPID HEMOGLOBIN & HEMATOCRIT RESULTS

	SIZE	REF#
AimStrip® Hemoglobin Starter Kit 200	200 tests	78200
<ul style="list-style-type: none"> ⇒ 1 AimStrip® Hemoglobin Meter ⇒ 200 Test Cartridges ⇒ 2 Optical Verifiers ⇒ 1 Carrying Case ⇒ User's Manual ⇒ 240 Capillary Tubes ⇒ 4 Code Chips ⇒ 4 AAA Batteries ⇒ 4 Results Recording Pads 		



STORE AT 2° to 30°C



AimStrip® Hemoglobin Meter	SIZE	REF#
	1 meter	78101

- ⇒ Includes: 1 Meter, 2 Optical Verifiers, 1 Carrying Case, 4 AAA Batteries



AimStrip® Hemoglobin Test Combo	SIZE	REF#
	50/60 combo	78250

- ⇒ 50 Test Cartridges
- ⇒ 60 Capillary Tubes
- ⇒ CLIA Waived



STORE AT 2° to 8°C



AimStrip® Hb Control Set	SIZE	REF#
	3 x 2mL	78103

- ⇒ 1 x 2mL Control Solution 0 (low)
- ⇒ 1 x 2mL Control Solution 1 (normal)
- ⇒ 1 x 2mL Control Solution 2 (high)



AimStrip® Hb Optical Verifiers	SIZE	REF#
	2 pack	78102

- ⇒ 2 Optical Verifiers



AimStrip® Capillary Transfer Tubes	SIZE	REF#
	50 pack	78155

- ⇒ 50 Capillary Tubes



AimStrip® Hb AC Adaptor	SIZE	REF#
	1 each	78111

- ⇒ 1 AC Adaptor (for meter)

CLIA
WAIVED



REF# 77301

AimStrip® TANDEM MEASURING SYSTEM

LIPID PROFILE & GLUCOSE MEASURING SYSTEM

METER USED FOR TESTING:

- AimStrip® Total Cholesterol Test Strips
- AimStrip® Triglycerides Test Strips
- AimStrip® Blood Glucose Test Strips

AimStrip® TANDEM BRAND
TEST STRIPS & CONTROL
SOLUTIONS ARE REQUIRED
BUT NOT INCLUDED IN KIT.



STORE AT 2° to 30°C



MULTIPLE PATIENT USE



1 KIT

INCLUDES:

- ➔ AimStrip® Tandem Meter
- ➔ Operation Manual
- ➔ Warranty Card
- ➔ Carrying Case
- ➔ 2 AAA Batteries
- ➔ Log Book




AimStrip® Tandem TEST STRIPS - For use with the AimStrip® Tandem Measuring System



**AimStrip® Tandem
TOTAL CHOLESTEROL TEST STRIPS**

SIZE	REF#
10 Strips	77310

- ➔ 10 Test Strips
- ➔ 10 Capillary Rods



**AimStrip® Tandem
TRIGLYCERIDES TEST STRIPS**

SIZE	REF#
10 Strips	77410

- ➔ 10 Test Strips
- ➔ 10 Capillary Rods



**AimStrip® Tandem
BLOOD GLUCOSE TEST STRIPS**

SIZE	REF#
50 Strips	77550

- ➔ Packed in 2 Bottles
- ➔ 25 Test Strips per Bottle

AimStrip® Tandem CONTROL SOLUTIONS

CONTROLS ARE FOR USE WITH AimStrip® TANDEM METER ONLY.

**AimStrip® Tandem TOTAL
CHOLESTEROL Control Solution**

SIZE	REF#
2 x 3mL	77312

- ➔ Control Solution Level 1 (3mL)
- ➔ Control Solution Level 2 (3mL)

**AimStrip® Tandem
TRIGLYCERIDES Control Solution**

SIZE	REF#
2 x 3mL	77412

- ➔ Control Solution Level 1 (3mL)
- ➔ Control Solution Level 2 (3mL)

**AimStrip® Tandem
GLUCOSE Control Solution**

SIZE	REF#
3 x 3.5mL	77512

- ➔ Control Solution Level 1 (3.5mL)
- ➔ Control Solution Level 2 (3.5mL)
- ➔ Control Solution Level 3 (3.5mL)



STORE AT 8° to 30°C

CLIA
WAIVED

OCCULT BLOOD FOB/GUAIAC

GUAIAC METHOD for detecting occult blood in stool



STORE AT 2° to 30°C



COMPLIANCE MONITORING
HIPAA Compliant program that reminds patients to return Collection Kits.

COMPLIANCE VERSION:

- ⇒ Full compliance version includes Compliance Monitoring
- ⇒ Return Labels and "Postage Pre-Paid" on Return Mailers
- ⇒ HIPAA compliant Patient Agreement forms

EACH KIT PROVIDES THE FOLLOWING:

- ⇒ On-Slide controls, detects Hb at 0.6mg of stool
- ⇒ Enhanced Readability results in 30 seconds
- ⇒ Collection Sticks & Developing Solution



FULL COMPLIANCE VERSION

The following kit **INCLUDES Compliance Monitoring and Pre-Paid Postage.**

Compliance Gold® FOB Test-ER, TRI-SLIDES

REF#
74252

ECONOMY VERSION

The following kits **DO NOT include Compliance Monitoring and Postage.**

Compliance Gold® FOB Test-ER, TRI-SLIDES

REF#
74250


Compliance Gold® FOB Test-ER, SINGLES

REF#
74200

DEVELOPER

	SIZE	REF#
Compliance Gold® FOB Developer	8 x 10 mL	74208

TRI-SLIDES:
50 TESTS PER KIT



SINGLES:
100 TESTS PER KIT



ENGLISH
AND
SPANISH
INSTRUCTIONS
AVAILABLE



**CLIA
WAIVED**

OCCULT BLOOD iFOB



STORE AT 2° to 30°C

30 TESTS PER KIT



IMMUNOLOGICAL TEST for occult blood in stool specimens



COMPLIANCE MONITORING
HIPAA Compliant program that reminds patients to return Collection Kits.

COMPLIANCE VERSION:

- ⇒ Full compliance version includes Compliance Monitoring
- ⇒ Return Labels & "Postage Pre-Paid" on Return Mailers
- ⇒ HIPAA compliant Patient Agreement forms



FULL COMPLIANCE VERSION

The following kit **INCLUDES Compliance Monitoring and Pre-Paid Postage.**

Compliance Gold® iFOB REF# 74034

- ⇒ Includes Test Devices, Collection Bottles & Collection Papers

ECONOMY VERSION

The following kit **DOES NOT include Compliance Monitoring and Postage.** (Return Mailers without postage are included).

Compliance Gold® iFOB Simple REF# 74030

- ⇒ Includes Test Devices, Collection Bottles & Collection Papers



STORE AT 2° to 30°C



- 15 Sets:
- ⇒ Return Mailers
 - ⇒ Collection Bottles
 - ⇒ Collection Papers

Patient Collection Kits

With Pre-Paid Postage REF# 74014

Use with Compliance Gold® iFOB

Without Postage REF# 74015

Use with Compliance Gold® iFOB Simple



Compliance Gold® iFOB Easy

SIZE REF#
30 tests 74033

- ⇒ Test Devices
- ⇒ Collection Bottles
- ⇒ No Compliance, Mailers or Postage



STORE AT 2° to 8°C



iFOB Controls

SIZE REF#
2 x 1mL 74211

- ⇒ Controls for occult blood tests
- ⇒ 1 x 1mL (Positive Control)
- ⇒ 1 x 1mL (Negative Control)

PREGNANCY TESTS



**CLIA
WAIVED
FOR
URINE**



STORE AT 2° to 30°C



AimStep® Pregnancy

SIZE	REF#
30 Tests	97030

- ⇒ Urine cassette method
- ⇒ Sensitive to 20 mIU/mL hCG
- ⇒ CLIA waived, also available in OTC



AimStep® Combo Pregnancy

SIZE	REF#
30 Tests	97730

- ⇒ Serum / Urine cassette method
- ⇒ Sensitive to 20 mIU/mL hCG
- ⇒ CLIA waived for urine testing



AimStrip® Pregnancy

SIZE	REF#
30 Tests	88530

- ⇒ Urine dipstick method
- ⇒ Sensitive to 20 mIU/mL hCG
- ⇒ CLIA waived, also available in OTC



STORE AT 2° to 30°C

AccuHome® PREGNANCY

OTC URINE MIDSTREAM METHOD



SIZE	REF#
1 Test	89101
2 Tests	89102

- ⇒ OTC final box packing
- ⇒ Sensitive to 25 mIU/mL hCG
- ⇒ CLIA waived



STORE AT 2° to 8°C

HCG CONTROLS

FOR USE WITH URINE PREGNANCY TESTS



SIZE	REF#
2 x 5mL	71002

- ⇒ Negative Control (1 x 5mL)
- ⇒ Positive Control (1 x 5mL)

FOR USE WITH SERUM PREGNANCY TESTS



SIZE	REF#
3 x 2.5mL	71003

- ⇒ Negative, Low Positive, and High Positive Control

**CLIA
WAIVED**

OVULATION TESTS



AimStep® Ovulation

SIZE	REF#
30 Tests	75030

- ⇒ Urine cassette method
- ⇒ OTC or professional use



AccuHome® Ovulation

SIZE	REF#
1 Test	75501

- ⇒ OTC Urine midstream method
- ⇒ Bulk packed

DIPSTICK METHOD FOR RAPID
DETECTION OF GROUP A STREP

StrepAim® STREP A



**CLIA
WAIVED**



- KITS INCLUDE:
- Dipsticks
 - Controls
 - Extract Reagents
 - Swabs
 - Tubes
 - Workstation

TEST DIRECTLY FROM
THROAT SWABS



STORE AT 2° to 30°C



StrepAim®

SIZE
25 Tests

REF#
73025



STORE AT 2° to 30°C



SERUM, PLASMA AND WHOLE BLOOD CASSETTE METHOD

AimStep® H. PYLORI

AimStep® MONO

**CLIA
WAIVED
FOR WHOLE
BLOOD**



SIZE
20 tests

REF#
33070



SIZE
25 tests

REF#
77025

KITS INCLUDE:



➤ H.Pylori Cassettes, Positive and Negative Controls, Developer Solution & Fingersticks

KITS INCLUDE:



➤ Mono Cassettes, Positive and Negative Controls, Developer Solution & Fingersticks



STORE AT 2° to 8°C

SEROLOGY



QUALITATIVE & SEMI-QUANTITATIVE PROCEDURES

Latex agglutination SLIDE METHOD

KITS INCLUDE: Latex, Controls, Buffer, Test Slide, and Pipettes

Aim™ ASO Test SIZE REF#
50 tests 31050

➤ for non-diluted serum to detect for Anti-streptolysin O

Aim™ CRP Test 100 tests 32100

➤ for non-diluted serum to detect for C-Reactive Protein

Aim™ RF Test 100 tests 35100

➤ for non-diluted serum to detect for Rheumatoid Factor

Carbon flocculation SLIDE METHOD

KIT INCLUDES: Carbon Antigen, Controls, Disposable Slides, and Pipettes

Aim™ RPR Test SIZE REF#
500 tests 36500

➤ for serum or plasma to detect for reagin antibodies

TABLET TESTS



STORE AT 2° to 30°C

AimTab™ BILIRUBIN TABLETS

REF# 13200
SIZE 100 tests

Tablet method for Bilirubin in urine.

- ⇒ Compare to Ictotest®, trademark by Siemens Healthcare
- ⇒ Test Mats included
- ⇒ Controls, see UA Liquid Controls Plus, Ref #51116



STORE AT 10° to 35°C

AimTab™ BLUING TABLETS

REF# 40021
SIZE 100 tests



Tablet method designed for use in all drug testing collections to prevent adulteration.

- ⇒ Meets Federal DOT Regulation 49 CFR 40.41
- ⇒ Not classified as an In-vitro Diagnostic test



STORE AT 2° to 30°C



AimTab™ KETONE TABLETS

Tablet method for the semi-quantitative determination of ketones (acetoacetic acid and acetone) in urine, serum, plasma & whole blood.

- ⇒ Detects as little as 5 mg/dL of acetoacetic acid in urine
- ⇒ Compare to Acetest®, trademark by Bayer Healthcare LLC.
- ⇒ Urine Control, see UA Liquid Controls Plus, Ref #51116
- ⇒ CLIA moderately complex (testing serum, plasma or whole blood)



Ketone Serum Controls

SIZE	REF#
2 x 2mL	13112

- ⇒ 1 Negative Control
- ⇒ 1 Positive Control
- ⇒ **Store at 2° to 8°C**

REF# 13100
SIZE 100 tests

AimTab™ REDUCING SUBSTANCES TABLETS



STORE AT 2° to 30°C

Tablet method for the semi-quantitative determination of the amount of reducing substances (glucose, lactose, fructose, galactose, pentoses) in urine, stool, or wine.

- ⇒ Wine making, used to estimate the residual sugar level of wine
- ⇒ Compare to Clinitest®, trademark by Bayer Healthcare LLC.
- ⇒ For Forensic Use Only



	SIZE	REF#
BOTTLED TABLETS	36 tests	13336
POUCHED TABLETS	20 tests	13320


- ⇒ 1 tablet per sealed foil pouch

TABLETING CONTRACT MANUFACTURING

A small-scale, personal contract tableting manufacturer. Focused on full development support for customers with small-scale manufacturing needs. cGMP compliant contract manufacturing services are available for any formulation. Compliant with Title 21 CFR part 11. CONTACT CUSTOMER SERVICE.



AimStrip® URINALYSIS STRIPS



STORE AT 2° to 30°C

PACKAGED IN BOTTLES
OF 100 REAGENT STRIPS

URINE DIPSTICKS

	REF#		REF#
AimStrip® US-G	50100	AimStrip® US-5	55100
⇒ for Glucose		⇒ for Glucose, Protein, Ketone, pH, Blood	
AimStrip® US-K	50121	AimStrip® US-5B	55200
⇒ for Ketone		⇒ for Blood, Glucose, Protein, Nitrite, Leukocytes	
AimStrip® US-2	52100	AimStrip® US-9	59100
⇒ for Glucose & Protein		⇒ for Glucose, Bilirubin, Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite	
AimStrip® US-3	53100	AimStrip® 10-SG	51100
⇒ for Glucose, Protein, pH		⇒ for Glucose, Bilirubin, Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite, Leukocytes	
AimStrip® US-4OB	54100		
⇒ for Glucose, Protein, Nitrite, Leukocytes			

PLEASE INQUIRE ABOUT OTHER AVAILABLE COMBINATIONS

UTI COMBO STRIPS



SIZE REF#
1 Test 52301

- ⇒ for Leukocytes and Nitrites
- ⇒ For OTC use

MICROALBUMIN STRIPS



AimStrip® MACRE Strips

SIZE REF#
25 Tests 52275

- ⇒ Sensitivity as low as 2.0 mg/dL
- ⇒ Detects Microalbumin & Creatinine as low as 10mg/L
- ⇒ Semi-Quantitative A:C ratio results in one minute
- ⇒ Same procedure as leading national brand



EXCLUSIVELY
READ ON
AimStrip® URINE
ANALYZER 2



MicroAlbumin 2-1 Combo Strips

25 Tests 52255

- ⇒ Urine Dipstick for Microalbumin and Creatinine
- ⇒ Determines Albumin to Creatinine ratio
- ⇒ For visual read only, analyzer not required

FOR THE SEMI-QUANTITATIVE DETECTION
OF ANALYTES IN URINE



AimStrip® URINE ANALYZER 2

REF# 51222



STORE AT 2° to 30°C



1 UNIT

- ⇒ Reads all combinations of strips
- ⇒ Single Test Mode, reads 1 strip at a time
- ⇒ Continuous Test Mode, reads 120 strips/hour
- ⇒ Automatic calibration for accurate results
- ⇒ Automatically saves the last 2000 results
- ⇒ Includes RS232C and USB ports to export results

LANGUAGES INCLUDED: ENGLISH, FRENCH, PORTUGUESE, SPANISH, ITALIAN



REF# 51110
Barcode reader included

AimStrip® URINE AUTO ANALYZER

LANGUAGES INCLUDED: ENGLISH, SPANISH, PORTUGUESE

- ⇒ High throughput, up to 500 tests per hour
- ⇒ Large touch screen for easy navigation
- ⇒ Reads 4, 8, 10 parameter strips
- ⇒ Automatic calibration for accurate results
- ⇒ Automatically saves the last 2000 results
- ⇒ Includes RS232C and USB ports to export results



STORE AT 2° to 8°C

UA CONTROLS

Liquid



	SIZE	REF#
UA Liquid Controls	2 x 5mL	51115
⇒ Negative Control (1 x 5mL), Positive Control (1 x 5mL)		
UA Liquid Controls	6 x 10mL	51161
⇒ Negative Control (3 x 10mL), Positive Control (3 x 10mL)		
UA Liquid Controls Plus	2 x 5mL	51116
⇒ Negative Control (1 x 5mL), Positive Control (1 x 5mL)		

DipTube



	SIZE	REF#
UA DipTube Controls	2 x 12mL	51112
⇒ Negative Control (1 x 12mL), Positive Control (1 x 12mL)		
UA DipTube Controls	4 x 12 mL	51142
⇒ Negative Control (2 x 12mL), Positive Control (2 x 12mL)		



STORE AT 2° to 30°C

PAPER

4 ROLLS/PK



REF#

Printer Paper	51224
⇒ Thermal paper to print 200 results per roll	
Printer Peel & Stick Paper	51234
⇒ Peel & stick printer paper to print 100 results/roll	
⇒ Conveniently place results into patient records	

BARCODE READER



REF#

Barcode Reader For Analyzers	51320
⇒ Use with AimStrip® Urine Analyzer 2 & AimStrip® Urine Auto Analyzer to record patient ID	



FOR
PROFESSIONAL
USE
ONLY
FOR PETS



VQPET H BLOOD GLUCOSE TESTS & SUPPLIES

FRESH WHOLE BLOOD (CAPILLARY, VENOUS, ARTERIAL)
CAN BE USED FOR TESTING FOR VETERINARY CLINICS.



CALIBRATED FOR
DOGS & CATS



GLUCOSE
DEHYDROGENASE-
FLAVIN ADENINE
DINUCLEOTIDE
(GDH-FAD) ENZYME
USED



VQPET H GLUCOSE STARTER KIT

REF# VGM30-431



1 KIT

INCLUDES:

- ⇒ 1 Pet Blood Glucose Meter
- ⇒ 10 Test Strips
- ⇒ 10 Sterile Lancets
- ⇒ 1 Lancing Device
- ⇒ 2mL Control Solution
- ⇒ Carrying Case
- ⇒ User's Manual
- ⇒ Quick Reference Guide
- ⇒ Logbook

TEST STRIPS



VQPET
GLUCOSE TEST STRIPS

REF#
VGS01-421

- ⇒ 50 Strips/Box
- ⇒ (2) 25 Strips/Vial

CONTROL SOLUTION

VQPET
CONTROL SOLUTION



REF#
VGM30-371

- ⇒ 1 x 2mL

Support Alternative Site Testing for Pets • Only 0.5µL Blood Sample Required from Pet



COVID-19 TESTS



STORE AT 2° to 30°C

QUALITATIVE & RAPID DETECTION OF CORONAVIRUS INFECTION

Status™ COVID-19/FLU A&B (SWAB)

REF# 33225
SIZE 25 tests

- ⇒ 1 test + 1 specimen = 3 Results
- ⇒ Rapid Immunoassay
- ⇒ Flocked Nasopharyngeal Swab
- ⇒ COVID-19 Sensitivity: 93.9%
- ⇒ COVID-19 Specificity: 100%
- ⇒ FDA Emergency Use Authorized (EUA)



- KITS INCLUDE:
- ⇒ 25 Test Devices
 - ⇒ 25 Extraction Reagent in Capsules
 - ⇒ 25 Sterile Swabs
 - ⇒ 1 Positive Control Swab
 - ⇒ 1 Negative Control Swab
 - ⇒ 1 Instructions for Use
 - ⇒ 1 Quick Reference Instructions

Celltrion DiaTrust™ COVID-19 AG RAPID TEST (SWAB)

REF# CT-P60 D-2 02
SIZE 25 tests

- ⇒ Lateral Flow Immunoassay
- ⇒ Nasopharyngeal Swab (NPS)
- ⇒ Clinical Sensitivity: 93.33% (28/30)
- ⇒ Clinical Specificity: 99.03% (102/103)
- ⇒ FDA Emergency Use Authorized (EUA)



- KITS INCLUDE:
- ⇒ 25 Test Devices
 - ⇒ 25 Sterile Swabs
 - ⇒ 25 Extraction Buffer
 - ⇒ 25 Filter Caps
 - ⇒ 1 Positive Control Swab
 - ⇒ 1 Negative Control Swab
 - ⇒ 1 Quick Reference Instructions



FaStep ANTIBODY RAPID DIAGNOSTIC TEST

REF# COV-W23
SIZE 20 tests

- ⇒ Combined Sensitivity: 100%
- ⇒ Combined Specificity: 98.8%
- ⇒ FDA Emergency Use Authorized (EUA)



- KITS INCLUDE:
- ⇒ 20 Tests
 - ⇒ 20 Disposable Pipettes
 - ⇒ 20 Sterile Safety Lancets
 - ⇒ 20 Alcohol Prep Pads
 - ⇒ Buffer
 - ⇒ Package Insert







www.germainelabs.com

FOR FURTHER INFORMATION CONTACT CUSTOMER SERVICE



 800.854.8446
210.692.4192
 210.692.4198 Fax

Hours 8:30am to 5:00pm
Central Standard Time

Errors Products received in error
must be reported within 10 days

Damaged Products damaged in transit
must be filed within 10 days

P.O.D Requests for proof of delivery
must be made within 60 days

Returns All returns must be authorized.
Products not returnable include
expired or short dated products,
opened products, and products
improperly handled or stored.

Remit to: PO Box 47015
San Antonio, Texas 78265

Shipping Address: 11030 Wye Drive
San Antonio, Texas 78217 USA



THE TESTS YOU WANT. THE ACCURACY YOU NEED!