



MICROALBUMIN 2-1 COMBO STRIPS

MicroAlbumin 2-1 Combo Strips are for the semi-quantitative detection of MicroAlbumin and Creatinine in urine. This test is for in-vitro diagnostic use.

CLIA
WAIVED

25 TESTS
PER BOTTLE



MicroAlbumin 2-1 Combo Strips
REF #52255

**Visual MicroAlbumin
and Creatinine Ratio
Made Easy!**

MicroAlbumin

Detects MicroAlbumin in concentration levels as low as 10mg/L

Creatinine

Detects Creatinine in concentration levels as low as 100mg/L



LOOK AT THESE GREAT FEATURES!



Same test procedure as the leading brand!



Determines the Albumin to Creatinine ratio.



Helps in the early detection of kidney disease in patients with diabetes or hypertension.



Increased sensitivity of Albumin as recommended by the ADA.



Store at room temperature.



⇒ CAP / CLSI formatted procedures available via fax or email
- just add your name!

CPT CODE: MicroAlbumin: 82044 QWReimbursement \$6.44*
Creatinine: 82570 QW.....Reimbursement \$7.28*

**THE TESTS YOU WANT.
THE ACCURACY YOU NEED!**

DISTRIBUTED BY:



www.germainelabs.com



800.854.8446
210.692.4192



210.692.4198 Fax